
Maestro 2.10 Manual

How to Power-Up the System

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Specifying Nuance Hardware Settings

➤ Hardware > Edit Hardware Setting > Camera

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Imaging the Cri Business card with Qdat Spots

- ✓ Acquiring an Image Cube
- ✓ Unmixing the Image Cube

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Monitor: Nothing (x64) [Live Stream]

File Edit Hardware Tools Window Help

LoadCube Zoom Full Image Pan ToolReset

ABURD Spectra Measure Display

Fluorescence Low Light Focus

Bring And Region Of Interest
 Settings: Max Sample Size: 822
 2x2 [0.2" x 4.4"] Max Full

Wavelength And Exposure
 Wavelength (nm): 580 Exposure (ms): 280.00

Color Scale: 500 800

Autopose Mono Autopose Cube

Filter/Wavelength Selection
 Start Stop End Narrow
 50 10 70

Size (50 x 10 x 70)

Focus Focus: 1

Multi Fiber Support
 Add

Lighting
 Open Shutter Turn On Illuminator

Flat Fielding
 Flat Field Acquire Flat Images

Sample ID

Notes

Acquire Mono Acquire Cube
 Acquire ByCS

35-B Cabot Road Woburn, MA 01801
 Tel: 781.935.9099 Fax: 781.935.3388
 Toll Free: 1.888.372.1242
 website: www.cri-inc.com
 email: [redacted]

Shrink Grow


Thumbnail Measurements Log Cursor

Selected	ROI Number	Cube ID	Avg Signal (counts)	Avg Signal (scaled counts)	Avg Signal (x10 ⁻⁶ photons/cm ² /s)	Avg Signal (DB)	Total Signal (counts)	Total Signal (scaled counts)

Hardware Status Protocol Cursor (C,Y) Average Signal Temperature
 Ready [redacted] [361, 185] [2863.908] Stage (C) [redacted] Chamber (C) [redacted]


Power On
 Run Continue Acquire, Acquire & Unload Stop

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
Acquiring Spectral Images

Acquiring Images Using Saved Saved Protocol

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Camera Binning, Zoom and Region of Interest

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
Specifying the Current Wavelength and Exposure

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➤ Adjusting the Current Wavelength

✓ To manually set the current wavelength

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➤ Adjusting the Exposure Time(Autoexpose)

❖ To autoexpose at the current wavelength

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❖ To use custom wavelengths and exposures

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Specifying Custom Wavelength and Exposure settings

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
➤ **To save custom collection settings as a table to load again later**

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➤ **To load a wavelength and exposure table**

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Making a Filter/Wavelength Selection

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Maestro Filter Setting Name	Excitation Filter (Range)	Emission Filter	Acquisition Settings*
White	Neutral attenuator	Transparent glass filter	500 to 800 in 10 nm steps
Blue	455 nm (435 to 480 nm)	490 nm longpass	500 to 720 in 10 nm steps
Green	523 nm (503 to 548 nm)	560 nm longpass	560 to 750 in 10 nm steps
Yellow	595 nm (576 to 621 nm)	635 nm longpass	630 to 800 in 10 nm steps
Orange	605 nm (586 to 631 nm)	645 nm longpass	640 to 820 in 10 nm steps
Red	635 nm (616 to 661 nm)	675 nm longpass	670 to 900 in 10 nm steps
Deep Red	681 nm (641 to 688 nm)	700 nm longpass	700 to 950 in 10 nm steps
NIR	704 nm (684 to 729 nm)	745 nm longpass	740 to 950 in 10 nm steps
ICG	735 nm (700 to 770 nm)	790 nm longpass	780 to 950 in 10 nm steps

* The ending wavelength can always be made longer if there is fluorescence beyond the suggested end point.

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
Acquisition Setup Using Multiple Filters

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
➤ Maintaining Compatible Exposure Time Ratio

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
➤ **Saving Many Filter Protocols in one Maestro Protocol**

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➤ **Taking a Reference Image for Flat Fielding**

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1. **Flat Field Check box**을 선택-Acquire Ref. Image은

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
Acquiring a Fluorescence Image Cube

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


- Low light Image(lower left)
- Low light Pose(second from left)
- Low light Composite(upper right)

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
Saving Low-Light Images

- **Saving the Composite**
File>Save Composite
- **Saving Images as Unscaled DATA**
File>Save All Images(As unscaled Data)

Assigning Sample IDs and Notes

➤ To Assign a Sample ID to an Image Cube

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
✓ To save Notes with an Image Cube

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Acquiring Timed Sequences of Image Cubes

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1. **Tools>Acquire Time Sequence**- Acquire Time Sequence dialog box 을 open 함.

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

- Save(As Displayed)
- Save All(As Displayed)
- Save image(As Un-scaled Data)
- Save image(As Un-scaled Data)

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
Viewing Cube Information

The screenshot shows a window titled 'Cube Info' with a scrollable area containing two tables. The first table lists general scan parameters, and the second table shows a grid of wavelength and exposure time data.

Cube Info	
Width (Pixels)	656
Height (Pixels)	520
Planes	23
Sample Id	
Comments	
Timestamp	2008-Jan-22 11:01:26.796976
mm/pixel	0.1510812
F Stop	4
Stage Height	1
Camera Gain	1
Camera Offset	64

Wavelength (nm)	500	510	520	530	540	550	560	570	580	590	600	610	620	6
Exposure Time (ms)	221.244	221.244	221.244	221.244	221.244	221.244	221.244	221.244	221.244	221.244	221.244	221.244	221.244	2

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Extracting an Image Plane from a Cube

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Unmixing Spectral Images

- **Opening a Spectral Library**
- **Opening an Image Cube**
- **Computing and Unmixing Spectra Automatically**
- **Computing and Unmixing Spectra Manually**
- **Saving An Unmixed Result Set**
- **Working With A Saved Result sets**
- **Checking Your Spectral Library**
- **Subtracting Spectra from a Cube**
- **Using Line Profiles to Analyze Signals**
- **Comparing Multiple Images**
- **Processing a Batch of Cubes**

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Computing and Unmixing Spectra Automatically

➤ **To compute Spectra using real component analysis:**

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Computing and Unmixing Spectra Manually

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✓**To hide a signal (e.g., background) in an unmixed image**

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✓**To restore default or saved colors**

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
✓**To add notes to the current spectral library**

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The screenshot displays the 'Manual Compute Spectra' window with three main sections:


- Known Spectrum (e.g. autofluorescence):** Shows a black curve on a graph with a y-axis from 0.0 to 1.0 and an x-axis from 580 to 700. A legend on the left shows C1 (red), C2 (green), and C3 (black).
- Mixed Spectrum (e.g. autofluorescence plus label):** Shows a red curve on a similar graph. The legend shows C1 (red), C2 (green), and C3 (black).
- Computed Spectrum:** Shows a red curve (Mixed Spectrum) and a blue dashed curve (Computed Spectrum) on a graph with a y-axis from 0 to 60 and an x-axis from 580 to 700. A 'Scaling' dropdown is set to 'Scale To Max'. Below the graph is a 'Fit Adjustment' section with 'Scale' checked, 'Fit Offset' unchecked, and an 'Error Scaler' slider.

At the bottom, there is a 'Transfer To Spectrum Number' dropdown set to '3', a 'Name' field containing 'C3', and a 'Color' dropdown set to 'Blue'. 'Transfer To Library' and 'Close' buttons are also visible.

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Saving An Unmixed Result Set

➤ To save the Current Unmixed result set

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

Using Line Profiles to Analyze Signals

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
Line Profiles drawn on component images

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Comparing Multiple Images

The screenshot displays the MicroManager software interface, specifically the 'Compare Images' window. The interface is divided into several panels:

- Left Panel:** Contains a 'Library' section with a 'Scaling' graph showing two curves (red and black) on a plot of intensity vs. wavelength (nm). Below the graph is a 'Notes' field and a color selection table with 10 entries (Red, Black, Blue, Yellow, Cyan, Magenta, Brown, White, Black, White).
- Top-Left Panel:** 'Grid Dimensions' section with fields for Rows (4), Cell Width (380), and Cell Height (283). It also includes 'Display Mapping' options for Gamma, Brightness, Contrast, Min/Max Dips, and Invert colors.
- Top-Right Panel:** 'Color Map' section with a 'Gray' color map and 'Overlays' section with 'Substrate Mask' set to Red.
- Center Panel:** A 2x2 grid of image windows, each titled 'Diss-1-4_0000'. The top-left window shows a dark image with two bright spots. The top-right window shows a similar image with a different processing. The bottom-left window shows the same image with a different processing. The bottom-right window shows the same image with a different processing. A red arrow points from the 'Diss-1-4_0000' window in the bottom-left of the grid to the 'Diss-1-4_0000' window in the bottom-left of the 'Thumbnail' section.
- Bottom Panel:** 'Thumbnail' section with tabs for 'Thumbnail', 'Measurements', 'Log', and 'Camera'. It shows five thumbnails: 'Diss-1-4', 'Diss-1-4_0000' (highlighted with a red border), 'Eye', 'AutoFL', and 'Camera'. Below the thumbnails are fields for 'Hardware Status' (Ready), 'Focal' (set to a value), 'Cursor (X,Y)', 'Average signal' (scaled counts/s), 'Temperature Stage (C)' (set to 30.0), and 'Temperature Chamber (C)' (set to 30.0). The 'Macros' section at the bottom right includes buttons for 'Run', 'Continue', 'Abort', 'Acquire & Upload', and 'Device...'

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
Quantifying Results

- **Measuring Regions**
- **Ignoring Smaller Regions**
- **Adjusting Region Transparency and Color**
- **Understanding Region Measurements**

Measuring Regions

➤ Finding Regions Automatically using Threshold Segmentation

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
➤ Manually Drawing and Modifying Regions

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
➤ **Dragging a Copy of a Region to Another Image**

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➤ **Copying all Regions to the Clipboard**

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
➤ Obtaining Accurate Measurement Data

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
Adjusting Region Transparency and Color

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Understanding Region Measurements

Selected	ROI Number	Cube ID	Avg Signal	Avg Signal/Exp.(ms)	Total Signal	Total Signal/Exp.(ms)	Max Signal	Area (Pbvals)	Area (mm ²)	Major Axis	Minor Axis
<input checked="" type="checkbox"/>	1		921.63	Unknown Exp.	427635.88	Unknown Exp.	1185.74	464.00	3.15	24.07	2.4
<input checked="" type="checkbox"/>	2		606.16	Unknown Exp.	248526.50	Unknown Exp.	1136.08	410.00	2.78	30.78	2.4
<input checked="" type="checkbox"/>	3		563.34	Unknown Exp.	171255.44	Unknown Exp.	1020.91	304.00	2.06	24.97	1.1

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
➤ Copying Measurement Data to the Clipboard

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➤ Saving measurement Data as a Text File

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
Customizing Spectral Displays

- **Adjusting Brightness and Contrast Levels**
- **Applying Overlays**
- **Adjusting a Cube's RGB Mapping**
- **Changing Components in composite Image**
- **Advanced Display Controls**
- **Displaying the Intensity Legend**

Adjusting Brightness and Contrast Levels

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Brightness ,Contrast, Red, Green
Blue

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
Changing Components in a Composite Image

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Advanced Display Controls

- ✓ **To open the Display Control utility**

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- ✓ **To link the Display control utility**

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
- ✓ To adjust a layer's brightness or contrast

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- ✓ To view images in fluorescence or brightfield

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
✓ **To insert commands into the new macro;**

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Co-localization Staining Detection

- Opening a Dataset for Co-localization Analysis
- Adjusting Threshold Mask Values
- Selecting Markers for Co-localization
- Saving and Loading Co-localization Settings
- Interpreting the Statistics
- Customizing the Statistics Display
- Drawing Regions of Interest
- Customizing the Composite Image Display
- Copying Images and Data
- Exporting Images and Data

b. The spectral components of the result set are listed here.

c. The thresholding values for each component can be adjusted here.

d. The markers and counterstains for co-localization can be selected here.

The screenshot shows the 'Co-localization' software interface. At the top, there's a 'Components' table with columns for Spectral Name, Marker, Color, and Visibility. Below it is a 'Thresholding' section with four rows, each for a marker (AF680, AF700, AF750, Autofluor...), including fields for Minimum Pixels, Peak Color, and a Threshold slider. The 'Co-localization' section at the bottom has checkboxes for 'Co-loc. Label', 'Show Co-loc.', 'Show All', 'Mask Color', and 'Visibility', along with two columns of checkboxes for 'Markers for Co-localization' and 'Denominator(Counterstain)'. On the right side, there's a 'Zoom' panel with buttons for 'Full Image', 'Fit', 'Red ROI', 'Area ROI', 'Cross ROI', and 'Clear ROIs'. A large image of a mouse is displayed in the center-right, with several colored spots (red, green, blue, yellow, white) overlaid on its body. A text box 'a. The 'Unmixed_Composite' image from the result set displays here.' points to this image. At the bottom right, there are buttons for 'Copy Data', 'Export Data', 'Copy Image', 'Export Image', and 'Close'. A text box 'e. Co-localization percentages for the selected markers and counterstains will display here.' points to the 'Statistics' area, which is currently empty. Another text box 'f. Use these buttons to copy and export data and images.' points to the bottom buttons.

Spectral Name	Marker	Color	Visibility
AF680	AF680	Red	<input type="checkbox"/>
AF700	AF700	Green	<input type="checkbox"/>
AF750	AF750	Blue	<input type="checkbox"/>
Autofluor...	Autofluor...	White	<input type="checkbox"/>

Marker Name	Minimum Pixels	Peak Color	Visibility
AF680	20	Red	<input type="checkbox"/>
Threshold			0.8000
AF700	20	Green	<input type="checkbox"/>
Threshold			0.8000
AF750	20	Blue	<input type="checkbox"/>
Threshold			0.8000
Autofluor...	20	White	<input type="checkbox"/>
Threshold			0.8000

Co-loc. Label	Show Co-loc.	Show All	Mask Color	Visibility
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Markers for Co-localization	Denominator(Counterstain)
<input type="checkbox"/> AF680	<input type="checkbox"/> AF680
<input type="checkbox"/> AF700	<input type="checkbox"/> AF700
<input type="checkbox"/> AF750	<input type="checkbox"/> AF750
<input type="checkbox"/> Autofluor...	<input type="checkbox"/> Autofluor...
	<input checked="" type="checkbox"/> All Image Pixels

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Selecting Markers for Co-localization

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The screenshot shows the 'Co-localization' software interface. The 'Components' section lists markers: AP500 (Red), AP700 (Green), AP750 (Blue), and Autofluor... (White). The 'Thresholding' section shows sliders for each marker's threshold. The 'Co-localization' section shows 'Overlap' as the label and 'yellow' as the mask color. The 'Statistics' section contains two tables.

	% Co-loc.	% Positivity AP500	% Positivity AP700	% Positivity AP750
PERCENT (%)	3.11%	9.66%	6.80%	6.50%

	% Positivity AP500	% Positivity AP700	% Positivity AP750
Within AP500	—	61.07%	24.93%
Within AP700	69.50%	—	28.52%
Within AP750	44.40%	50.16%	—

Saving and Loading Co-localization Settings

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✓ To Save a Co-localization setting file

- 1.co-localization window을 열어둔 채, **File>Save Co-localization**
2. save 버튼을 클릭함.

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3. 선택된 settings file은 Setting box에 display됨.

✓ **Co-localization settings must be compatible:**

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	% Positivity: AF680	% Positivity: AF700	% Positivity: AF750	Pixel Count
Within AF680	e --	61.07%	34.40%	14391
Within AF700	66.60%	--	38.82%	13195
Within AF750	48.49%	50.18%	--	10208

	ROI Number	Marker Area (pixels)	Total Signal (counts)	Total Signal (scaled counts/μ)	Avg Signal (counts)	Avg Signal (scaled counts/μ)	Max Signal (counts)	Max Signal (scaled counts/μ)
AF680	f Overlap	4638	340378.19	29.10	73.39	0.00628	139.53	0.0119
AF680	Full Image	14391	1579186.88	135.03	4.36	0.000373	180.32	0.0154
AF700	Overlap	4638	375131.28	32.08	80.88	0.00692	146.45	0.0125
AF700	Full Image	13195	1337325.88	114.35	3.70	0.000316	162.67	0.0139
AF750	Overlap	4638	208150.08	17.80	44.88	0.00384	74.00	0.00624
AF750	Full Image	10208	943201.06	80.65	2.61	0.000223	314.24	0.0269
AutoFluor...	Overlap	4638	957427.94	81.86	206.43	0.0177	424.75	0.0363
AutoFluor...	Full Image	148922	22637338.00	1935.59	62.55	0.00535	775.07	0.0663

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Customizing the Statistics Display

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The screenshot shows a software window titled 'Co-localization' with a 'Statistics' tab selected. The main window displays a grayscale image of a mouse with several colored regions (red, green, blue) overlaid. The left sidebar contains settings for components (AF680, AF750, AF750, Autofluor...) and checkboxes for displaying Co-localization, Peak Intensity/Area, In-Channel, and Component Statistics. The bottom right corner features buttons for 'Copy Data', 'Export Data', 'Copy Image', 'Export Image', and 'Close'.

Statistics Panel Data:

	% Co-loc.	% Positivity: AF680	% Positivity: AF750	% Positivity: AF750	Denominator
Percent (%)	3.11%	9.66%	6.66%	6.67%	
Pixel Count	1528	13191	13195	10306	148922

	% Positivity: AF680	% Positivity: AF750	% Positivity: AF750	Pixel Count
Within AF680	—	61.07%	34.40%	34391
Within AF750	66.60%	—	30.82%	13195
Within AF750	48.49%	50.10%	—	10306

	RDE Number	Marker Area (pixels)	Total Signal (counts)	Total Signal (scaled counts/a)	Avg Signal (counts)	Avg Signal (scaled counts/a)	Max Signal (counts)	Max Signal (scaled counts/a)
AF680	Overlap	4536	340070.19	25.10	73.39	0.00628	139.53	0.0119
AF680	I	14391	1578186.88	105.05	43.36	0.000773	180.32	0.0154
AF750	Overlap	4536	375131.29	22.08	80.88	0.00092	146.45	0.0125
AF750	I	13195	1337325.88	114.58	2.70	0.00016	162.87	0.0139
AF750	Overlap	4536	208159.98	17.88	44.88	0.00384	74.09	0.00634
AF750	I	10306	942001.05	80.65	2.61	0.00023	214.34	0.0269
Autofluor...	Overlap	4536	957427.04	61.86	206.43	0.0177	424.75	0.0363
Autofluor...	I	148922	22837338.00	1538.59	62.55	0.00335	775.07	0.0663

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➤ Applying Regions from a Saved ROI file

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Customizing the Composite Image Display

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

Dynamic Contrast Enhancement(Dyce)

- **How Dyce Works**
- **Preparing for Dyce imaging**
- **Acquiring a Monochrome Dyce Dataset**
- **Acquiring a Multispectral Dyce Dataset**
- **Recalling/Opening a Saved Dyce Cube Sequence**
- **Creating a Dyce Time Series**
- **Displaying/Playing Dyce Dataset**


How Dyce Works

1. Inject a bolus of NIR dye*Indocyanine Green(ICG)을 권장함
2. Dynamic series of image 사진을 찍음
3. Unmix

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
➤ Selecting Binning and ROI

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➤ Selecting a Collection Type

- **Mono(broadband emission)**-monochrome image을 정해진 time interval 로 acquire함
- **Mono(using current Acquire wavelength)**
- **Cube(using Current Acquire settings)**

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Collection Type

Mono (broadband emission)
 Mono (using current Acquire wavelength)
 Cube (using current Acquire settings)

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
Recalling/Opening a Saved Dyce Cube Sequence

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4. Create Dyce Time Series From A sequence of Cube 버튼을 클릭함.
 - Reloading a Time Series' Parent Cube Sequence
 1. Dyce time series cube을 열어놓채, Load Cube From This Sequence 버튼 클릭함
 2. Choose Cube in Sequence dialog에서 cube을 오픈함

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
➤ Changing the Display Mapping

➤ Selecting a Color Map

➤ Displaying Dimension-Reduced Data